



User Manual BatInspector

SoftwareVersion: **0.2**

Date: **23-05-30**

BatInspector is a software to analyze and manage recordings of bat calls.
It supports Elekon BatCorder projects as well as plain WAV files.

Main Features:

- manage bat recordings
- Analyze bat calls:
 - frequency max
 - frequency min
 - duration
 - many more
- Waterfall diagrams:
 - zoom in and out
 - zoom to detected call
- play recordings in original and 10x stretched speed
- Automatic species detection (experimental)
- Manual species detection
- Powerful data filtering
 - user defined filter expressions
- Generate reports



Index

1 Document History.....	2
2 Introduction.....	3
3 Installation and setup.....	3
3.1 Installation.....	3
4 Basic Workflow.....	4
4.1 Copy Data to Workstation.....	4
4.2 Open a Project.....	4
4.3 Analyze calls automatically.....	4
4.4 Working with Filters.....	4
4.5 Analyze calls manually.....	4
5 Menues and Screens.....	4
6 Scripting Language.....	5
6.1 Syntax.....	5
6.2 Commands.....	5
6.2.1 SET.....	5
6.2.2 CALL.....	5
6.2.3 FOR ... END.....	6
6.2.4 WHILE ... END.....	6
6.2.5 BREAK.....	7
6.2.6 LOG.....	7
7 Expression Parser Syntax.....	8
7.1 Operators.....	8
7.2 Data Types of Expressions.....	9
7.3 Supported Functions.....	9
7.3.1 Trigonometrical Functions.....	9
7.3.2 Mathematical Functions.....	10
7.3.3 String functions.....	11
7.3.4 Miscellaneous Functions.....	11
7.3.5 Functions to Handle CSV Files.....	12
7.3.6 Functions to deal with project informations.....	13
7.3.7 Functions to deal with WAV files.....	15

1 Document History

Date	Author	Changes
22-03-15	CMU	Initial version



2 Introduction

3 Installation and setup

3.1 Installation

- Install python 3.10
<https://www.python.org/downloads/>
It has to be installed in the location [C:\Program](#) Files\python310
- Install tensorflow 2
<https://www.tensorflow.org/install>
- Install BatInspector:
call setup.exe and follow the instructions




4 Basic Workflow

4.1 Copy Data from recorder to PC

In order to work with the recorded data, they must first be copied to the PC. Copy the data to a folder in your hard disk. Each project or recording session should have its own directory.

4.2 Open a Project



Select a directory in the file browser on the left side of the main window and click the directory containing the project to open.

Directories containing projects with project files are marked in cyan or green color.



4.2.1 Format of the project file

BatInspector is able to read project files in the Elekon data format.

The project file must have the extension *.bpr and the recordings must be located in a subfolder named „Records“:

	Organisieren	Neu
2022 > 20220719		
Name	Änderungsdatum	
log	27.08.2022 18:26	
Records	29.08.2022 21:05	
20220719.bpr	29.08.2022 21:06	
report.csv	29.08.2022 21:06	

The project file has an XML format, containing a list with all the files and some additional information:

```
1 <?xml version="1.0" encoding="utf-8"?>
2 <BatExplorerProjectFile xmlns:xsi="http://www.w3.org/2001/XMLSchema"
3   <Name>20220719</Name>
4   <Type>Elekon</Type>
5   <Records>
6     <Record File="20220719_0001.wav" Name="20220719_0001" />
7     <Record File="20220719_0002.wav" Name="20220719_0002" />
8     <Record File="20220719_0003.wav" Name="20220719_0003" />
9     <Record File="20220719_0004.wav" Name="20220719_0004" />
10    <Record File="20220719_0005.wav" Name="20220719_0005" />
```




4.2.2 Format of the Report File



4.3 Analyze calls automatically

4.4 Working with Filters



4.5 Analyze calls manually

5 Menues and Screens



6 Scripting Language

6.1 Syntax

The syntax of each script line is:

`<command> <Par1> <Par2> ... <Parn>`

Each line contains a command with parameters.

A parameter may be

- a fixed value: `22`
- a string enclosed in „“: `„this is a sample string“`
- an expression: `=(cos(a) + 5)`

An expression is surrounded by the start character „=(„ and the end character „)“:

`SET a =(cos(b) * 2)`

A detailed description of the expression syntax is provided in chapter 7

Before and between command and parameters arbitrary white space (space or tab) is allowed.

A line starting with the character „#“ is considered as a comment

`# this is a comment and will not be executed`

6.2 Commands

6.2.1 SET

Function: Assigns a value to a variable. If the variable does not exist, it will be created

Syntax: `SET <VarName> <Value>`

Parameter: VarName: Name of the variable
Value: Value to assign the variable

Example: `SET a 33`
`SET b =(a + 5)`

6.2.2 CALL

Function: calls a script. All variables of the calling script are accessible in the called script.

Syntax: `CALL <FileName>`

Parameter: FileName: Name of the script file (full path or only name. If only the name is provided, the file must be located in the directory <ScriptDir> from the application settings.

Example: `CALL C:\scripts\script1.scr`
`CALL script2.scr`



6.2.3 FOR ... END

Function: For Loop

Syntax: FOR <iteratorName> <start> <end>
 # a block of code
:
END

Parameter: iteratorName: name of the iterating variable
start: start value for the iterator
end: end value for the iterator

Example:

```
FOR i 1 10  
    LOG i  
END
```

6.2.4 WHILE ... END

Function: While Loop

Syntax: WHILE <expression>
 # a block of code
:
END

Parameter: expression: an expression with boolean result. As long the expression is true, the loop will be executed

Example:

```
SET i 0  
WHILE =(i < 10)  
    LOG i  
    SET i =(i + 1)  
END
```



6.2.5 BREAK

Function: breaks the execution of a loop (applicable for FOR and WHILE loops)

Syntax: BREAK

Parameter: none

Example:

```
WHILE =(a < 100)
  LOG a
  IF =(b >0) && ( a > 50)
    BREAK
  END
  SET a =(a + 1)
END
```

6.2.6 LOG

Function: write a message to the log window

Syntax: LOG <message> <type>

Parameter: message: the message to log (fixed val. Or expression)

type: optional parameter for message type:

INFO (default)

ERROR

WARNING

DEBUG

Example:

```
LOG „this should never happen“ ERROR
```



7 Expression Parser Syntax

With the integrated formula parser it is possible to express complex logic formulas or mathematical calculations.

This is useful to create filter expressions and scripts.

7.1 Operators

Operator	Description
+	Addition
-	Subtraction
*	Multiplication
/	Division
&&	Logical AND
	Logical OR
!	Logical NOT
<	Lesser than
<=	Lesser or equal than
>	Greater than
>=	Greater or equal than
==	Equal
!=	Not equal
()	Braces for math terms or function arguments
[]	Square braces for array index
„“	Double quotes to enclose strings



7.2 Data Types of Expressions

Data Type	Range	Examples
RT_BOOL	True, False	(a <= 3.0)
RT_INT64	$-2^{63} \dots 2^{63}$	-236, 224
RT_UINT64	$0 \dots 2^{64}$	23456
RT_HEX	0x0 ... 0xFFFFFFFFFFFFFFFF	0xFA10
RT_FLOAT	64 bit floating point number	-3.5, 66.8767, -3.23E-12
RT_COMPLEX	Complex numbers	1+2.0i, -4.5 – 3i, 2.3 + 1i
RT_STR	String	„this is a string“, „1.0“

7.3 Supported Functions

7.3.1 Trigonometrical Functions

sin(x)

Function returns sinus of a number",
 Params x: a number
 Returns sinus of x

asin(x)

Function returns arcus sinus of a number
 Params x: a number
 Returns arcus sinus of x

sinh(x)

Function returns sinus hyperbolicus of a number",
 Params x: a number
 Returns sinus hyperbolicus of x

cos(x)

Function returns cosinus of a number",
 Params x: a number
 Returns cosinus of x

acos(x)

Function returns arcus cosinus of a number",
 Params x: a number
 Returns arcus cosinus of x

cosh(x)

Function returns cosinus hyperbolicus of a number",
 Params x: a number



Returns cosinus hyperbolicus of x

tan(x)

Function returns tangens of a number",

Params x: a number

Returns tangens of x

atan(x)

Function returns arcus tangens of a number",

Params x: a number

Returns arcus tangens of x

7.3.2 Mathematical Functions

sqrt(x)

Function returns square root of a number

Params x: a number

Returns square root of x

pow(x, e)

Function calculate power of a number",

Params x: a number

e: exponent

Returns x^e (x to the power of e)

exp(x)

Function calculate power of e",

Params x: a number

Returns e^x (e to the power of x)

ln(x)

Function calculate the natural logarithm of a number

Params x: a number

Returns natural logarithm of x

abs(x)

Function calculate the absolute value of a number (real or complex)

Params x: a number

Returns absolute value of x

arg(x)

Function calculate the argument of a complex number

Params x: a complex number

Returns arg x



7.3.3 String functions

strlen(s)

Function calculate the length of a string

Params s: a string

Returns length of string s

strListCnt(str)

Function counts the number of items in a ,;‘ separated list of strings

Params str: a string

Returns number of items in list

strListItem(str, item)

Function returns the selected item in a ,;‘ separated list of strings

Params str: a string

 item: index of selected item (0..n)

Returns selected item of list

substr(str, start, length)

Function calculate the substring of a string

Params str: a string

 start: start character (0... <length of str – 1>)

 length: length of substr (0... <length of str – 1 - start>)

Returns the sub string

toupper(str)

Function convert string to upper case

Params str: a string

Returns the string in upper case



7.3.4 Miscellaneous Functions

cast(exp, type)

Function cast an expression to a specific type

Params exp: an expression

type: a type

(„RT_BOOL“, „RT_FLOAT“, „RT_INT64“, „RT_UIN64“, „RT_STR“, „RT_HEX“
„RT_COMPLEX“)

Returns length of string s

if(exp, argTrue, argFalse)

Function evaluate boolean expression and return corresponding argument

Params exp: a boolean expression

argTrue: an expression for exp == true

argFalse: an expression for exp == false

Returns argTrue or argFalse depending on result of exp

vars()

Function return a list of variables and their values

Params none

Returns list of all defined variables and their values

consts()

Function return a list of constants and their values

Params none

Returns list of all defined constants and their values



7.3.5 Functions to Handle CSV Files

openCsv(fileName,withHeader,separator)

Function: opens a csv file and reads it to memory

Parameter: fileName: Name of the script file (full path or only name. If only the name is provided, the file must be located in the directory <ScriptDir> from the application settings.name of the file (full path)
withHeader: (optional, default 0) 1: first row contains header with column names. The columns can be accessed via the names if opened in this mode.
Separator (optional, default=";") separation character in lines for the cells

Returns: >=0: a handle to access the file, <0: error opening the file

closeCsv(handle, write)

Function: closes a csv file and releases the memory"

Parameter: handle: file handle to access csv file
write: (optional, default:0) 1: write file to disk

Returns nothing

getCell(handle, row, col)

Function: get cell content of a previously opened file

Parameter: handle: file handle to access csv file
row row number (1..n)
col col name (if opened with option ,withHeader') or col number(1..n)

Returns content of cell as string

setCell(handle, row, col, value)

Function: set a cell content of a previously opened csv file

Parameter: handle: file handle to access csv file
row row number (1..n)
col col name (if opened with option ,withHeader') or col number(1..n)
value value to write to cell (string)

getRowCount(handle)

Function: get row count of a previously opened csv file (including header)

Parameter: handle: file handle to access csv file

Returns: number of rows in opened file



7.3.6 Functions to deal with projects

createPrjFromFiles(name, srcDir, destDir, maxFiles, maxFileLen, latitude, longitude, landscape, weather)

Function: create a project from a list of WAV files. If the number of WAV files exceeds the specified maximum number of files the projects gets split into multiple projects. The name of the subprojects are supplemented with an index nr: <prjName_<nnn>. Files that exceed the specified maximum file length are split into several file names. The resulting file names are supplemented with index numbers: <fileName>_<nnn>.wav

Parameter: name: name of the project (without extension)
srcDir: source folder containing WAV files
destDir: destination folder containing project folders
maxFiles: maximum number of files per project
maxFileLen: maximum file length per WAV file [s] (float)
latitude: geographical position
longitude: geographical position
landscape: short description of the landscape of the recording location
weather: short description of weather conditions during recording

Returns: error code: 0=OK

inspectPrj()

Function: identify calls, parameters of calls and the species in the recordings of the currently opened project. The results of the inspection process are stored in a report file named „report.csv“. The result file is stored in a sub folder with the name of the selected modell and is located in the project folder.

Parameter: none

Returns: error code: 0=OK

openPrj(directory, name)

Function: open a project.

Parameter: directory: directory that contains the project
name: project name (*.bpr)

Returns: projectList: a list of ,;‘ separated project names



7.3.7 Functions to deal with project informations

calcProbabilityRatios()

Function: calculate the ratio for the probabilities of classification for the 1st and 2nd rank for each call in the currently opened project.

Parameter: none

Returns: nothing

getCallCount(index)

Prerequisite: a record file must be generated first

Function: get number of calls for a specified file in currently opened project

Parameter: index: index of file (0..n)

Returns: number of calls

getCallInfo(fIndex, cIndex, infoType)

Prerequisite: a record file must be generated first

Function: get call information for a specified call in a specified file in currently opened project

Parameter: fIndex: index of file (0 ... n)

cIndex: index of call (0 ... m)

infoType: type of information

„SPEC_AUTO“: automatically detected species, classifier 1

„SPEC_MAN“: manually detected species

„PROB_RATIO“: ratio betw. 1st and 2nd rank of classification

„PROBABILITY“: probability for classified species

Returns: the requested information

getPrjFileCount()

Function: get nr of files in open project

Parameter: none

Returns: number of files in currently opened project

getFileName(index)

Function: get full path file name of file in project opened project

Parameter: index: index of file (0... n)

Returns: full path of wav file



getFileInfo(index, infoType)

Prerequisite: a record file must be generated first

Function: get file specific information in currently opened project

Parameter: index: index of file (0... n)
infoType: type of information to set
 „NAME“: name of the file
 „DURATION“: duration of the recording [s]
 „SAMPLE_RATE“: sampling rate [Hz]

Returns: the specified data

setFileInfo(index, infoType, data)

Prerequisite: a record file must be generated first

Function: set file specific information in currently opened project

Parameter: index: index of file (0... n)
infoType: type of information to set
 „SELECT“: set the state of the checkbox „Select“
 („TRUE“ or „FALSE“)

data: the data
Returns: error code

setCallInfo(fIndex, cIndex, infoType, data)

Prerequisite: a record file must be generated first

Function: set call information for a specified call in a specified file in currently opened project

Parameter: fIndex: index of file (0 ... n)
cIndex: index of call (0 ... m)
infoType: type of information
 „SPEC_MAN“: manually detected species
data: the data

Returns: the requested information

getNrOfSpecies(fIndex)

Prerequisite: a record file must be generated first

Function: get number of auto detected species for spec. file in open project

Parameter: fIndex: index of file (0 ... n)

Returns: nr of automatically detected species



getRankSpecies(fIndex, rank)

Prerequisite: a record file must be generated first

Function: get the species name for the specified rank in specified file in open project

Parameter: fIndex: index of file (0 ... n)
rank: rank (1 ... n, 1= most calls)

Returns: species name for specified rank

getRankCount(fIndex, rank)

Prerequisite: a record file must be generated first

Function: get the nr of calls of species for the specified rank in specified file in open project

Parameter: fIndex: index of file (0 ... n)
rank: rank (1 ... n, 1= most calls)

Returns: nr species in recording for specified rank



7.3.8 Functions to deal with WAV files

bandPass(file, fMin, fMax)

Function: sets the sampling rate of a file

Parameter: file: name of the file (full path) or alternatively file index of file in currently open project
fMin: low frequency limit for band pass [Hz]
fMax: high frequency limit for band pass Hz]

Returns: error code

getSampleRate(fileName)

Function: returns the sampling rate of a file

Parameter: fileName: name of the file (full path)

Returns: the sampling rate

rescaleSampleRate(fileName, factor)

Function: multiplies the sampling rate of a file with a factor

Parameter: fileName: name of the file (full path)
factor: factor to multiply sampling rate

Returns: the new sampling rate

SetSampleRate(fileName, sampleRate)

Function: sets the sampling rate of a file

Parameter: fileName: name of the file (full path)
sampleRate: new value for sampling rate

Returns: nothing